



COT Rapid Test Panel (Urine) Package Insert

REF DCT-114 English

A rapid test for the qualitative detection of Cotinine (nicotine metabolite) in human urine.
For determination of smoking status only. Not intended for medical diagnostic use.

INTENDED USE

The COT Rapid Test Panel (Urine) is a rapid chromatographic immunoassay for the detection of Cotinine in human urine at a cut-off concentration of 200ng/mL. This test will detect other related compounds, please refer to the Analytical Specificity table in this package insert.

This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography and mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

SUMMARY

Cotinine is the first-stage metabolite of nicotine, a toxic alkaloid that produces stimulation of the autonomic ganglia and central nervous system when in humans. Nicotine is a drug to which virtually every member of a tobacco-smoking society is exposed whether through direct contact or second-hand inhalation. In addition to tobacco, nicotine is also commercially available as the active ingredient in smoking replacement therapies such as nicotine gum, transdermal patches and nasal sprays.

In a 24-hour urine, approximately 5% of a nicotine dose is excreted as unchanged drug with 10% as cotinine and 35% as hydroxycotinine; the concentrations of other metabolites are believed to account for less than 5%. While cotinine is thought to be an inactive metabolite, it's elimination profile is more stable than that of nicotine which is largely urine pH dependent. As a result, cotinine is considered a good biological marker for determining nicotine use. The plasma half-life of nicotine is approximately 60 minutes following inhalation or parenteral administration.² Nicotine and cotinine are rapidly eliminated by the kidney; the window of detection for cotinine in urine at a cutoff level of 200 ng/mL is expected to be up to 2-3 days after nicotine use.

The COT Rapid Test Panel (Urine) is a rapid urine screening test that can be performed without the use of an instrument. The test utilizes a monoclonal antibody to selectively detect elevated levels of Cotinine in urine. The COT Rapid Test Panel (Urine) yields a positive result when the Cotinine in urine exceeds 200 ng/mL.

PRINCIPLE

The COT Rapid Test Panel (Urine) is an immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody.

During testing, a urine specimen migrates upward by capillary action. Cotinine, if present in the urine specimen below 200 ng/mL, will not saturate the binding sites of antibody coated particles in the test. The antibody coated particles will then be captured by immobilized Cotinine conjugate and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the Cotinine level exceeds 200 ng/mL because it will saturate all the binding sites of anti-Cotinine antibodies.

A drug-positive urine specimen will not generate a colored line in the test line region because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that that proper volume of specimen has been added and membrane wicking has occurred.

REAGENTS

The test contains mouse monoclonal anti-Cotinine antibody-coupled particles and Cotinine-protein conjugate. A goat antibody is employed in the control line system.

PRECAUTIONS

- For medical and other professional *in vitro* diagnostic use only. Do not use after the expiration date.
- The test should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test should be discarded according to local regulations.

STORAGE AND STABILITY

Store as packaged at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch or label of the closed canister. The test must remain in the sealed pouch or closed canister until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

SPECIMEN COLLECTION AND PREPARATION

Urine Assay

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed settle to obtain a clear supernatant for testing.

Specimen Storage

Urine specimens may be stored at 2-8°C for up to 48 hours prior to assay. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

MATERIALS

Materials Provided

- Test Panels
- Package insert


Materials Required But Not Provided

- Specimen collection container
- Timer

DIRECTIONS FOR USE

Allow the test, urine specimen, and/or controls to reach room temperature (15-30°C) prior to testing.

- Remove the test panel from the sealed pouch and use it within one hour.
- Remove the cap.



- With the arrow pointing toward the urine specimen, immerse the test panel vertically in the urine specimen for at least 10 to 15 seconds. **Immerse the strip to at least the level of the wavy lines, but do not touch the plastic device.**
- Replace the cap and place the test panel on a non-absorbent flat surface.
- Start the timer and wait for the colored line(s) to appear.
- The result should be read at 5 minutes. Results may be stable up to 1 hour after test initiation.

INTERPRETATION OF RESULTS

(Please refer to the illustration above)

NEGATIVE: Two lines appear. One colored line should be in the control line region (C), and another apparent colored line should be in the test line region (T). This negative result indicates that the Cotinine concentration is below the detectable level (200 ng/mL).

NOTE: The shade of color in the test line region (T) may vary, but it should be considered negative whenever there is even a faint colored line.

POSITIVE: One colored line appears in the control line region (C). No line appears in the test line region (T). This positive result indicates that the Cotinine concentration exceeds the detectable level (200 ng/mL).

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test. If the problem persists, discontinue using the lot immediately and contact your local distributor.

QUALITY CONTROL

A procedural control is included in the test. A colored line appearing in the control line region (C) is considered an internal positive procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

- The COT Rapid Test Panel (Urine) provides only a qualitative, preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.^{1,2}
- It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
- Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
- A positive result indicates presence of the drug or its metabolites but does not indicate level of intoxication, administration route or concentration in urine.
- A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
- Test does not distinguish between drugs of abuse and certain medications.

EXPECTED VALUES

This negative result indicates that the Cotinine concentration is below the detectable level of 200ng/mL. Positive result means the concentration of Cotinine is above the level of 200ng/mL. The COT Rapid Test Panel has a sensitivity of 200ng/mL.

PERFORMANCE CHARACTERISTICS

Accuracy

A comparison was conducted using the COT Rapid Test Panel (Urine) and GC/MS. The following results were tabulated:

Method	GC/MS		Total Results
	Positive	Negative	
COT Rapid Test Panel	88	4	92
	91	155	250
% Agreement	98.7%	97.5%	

Analytical Sensitivity

A drug-free urine pool was spiked with Cotinine at the following concentrations: 0 ng/mL, 100 ng/mL, 150 ng/mL, 200 ng/mL, 250 ng/mL, 300 ng/mL and 600 ng/mL. The results demonstrate > 99% accuracy at +50% above and 50% below the cut-off concentration. The data are summarized below:

Cotinine Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0	30	30	0
100	-50%	30	30	0
150	-25%	30	27	3
200	Cut-off	30	15	15
250	+25%	30	4	26
300	+50%	30	0	30
600	+300%	30	0	30

Analytical Specificity

The following table lists compounds that are positively detected in urine by the COT Rapid Test Panel (Urine) at 5 minutes.

Compound	Concentration (ng/mL)
(-)Cotinine	200
(-)Nicotine	5,000

Precision

A study was conducted at three hospitals by laypersons using three different lots of product to demonstrate the within run, between run and between operator precision. An identical panel of coded specimens containing, according to GC/MS, no Cotinine, 25% Cotinine above and below the cut-off, and 50% Cotinine above and below the 200 ng/mL cut-off was provided to each site. The results are given below:

Cotinine Concentration (ng/mL)	n per site	Site A		Site B		Site C	
		-	+	-	+	-	+
0	10	10	0	10	0	10	0
100	10	10	0	10	0	10	0
150	10	9	1	9	1	9	1
250	10	1	9	1	9	2	8
300	10	0	10	0	10	0	10

Effect of Urinary Specific Gravity

Fifteen urine specimens of normal, high, and low specific gravity ranges were spiked with 100 ng/mL and 300 ng/mL of Cotinine. The COT Rapid Test Panel (Urine) was tested in duplicate using the fifteen neat and spiked urine specimens. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

Effect of Urinary pH

The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with Cotinine to 100 ng/mL and 300 ng/mL. The spiked, pH-adjusted urine was tested with the COT Rapid Test Panel (Urine) in duplicate. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or Cotinine positive urine. The following compounds show no cross-reactivity when tested with the COT Rapid Test Panel (Urine) at a concentration of 100µg/mL.

Non Cross-Reacting Compounds

4-Acetamidophenol	4-Dimethylaminoantipyrine	Lithium carbonate	Phentermine
Acetone	Diphenhydramine	Loperamide	trans-2-Phenyl
Acetophenetidin	5,5-Diphenylhydantoin	Maprotiline	cyclopropylamine
Acetylsalicylic acid	Sopramide	Meprobamate	l-Phenylephrine
N-Acetylprocainamide	Doxylamine	Methamphetamine	β-Phenylethylamine
Albumin	Egonine	Meprobamate	Phenylpropranolamine
Aminopyrine	Egonine methylester	Methadone	(d,l-norephedrine)
Amtripyline	EDDP	d-Methamphetamine	(±) Phenylpropranolamine
Amobarbital	EFavirenz (Sustiva)	l-Methamphetamine	Prednisolone
Amoxapine	EMDP	Methaqualone	Prednisone
Amoxicillin	Epheдрine	Methoxyphenamine	5β-Pregnone-3α, 17α, 21-triol
l-Amphetamine	l-Ephedrine	(-) 3,4-Methylenedioxy-amphetamine (MDA)	Procaine
Ampicillin	(±)-Epinephrine	(+) 3,4-Methylenedioxy-methamphetamine	Promazine
Apomorphine	l-Epinephrine	(MDMA)	Promethazine
l-Ascorbic acid	Erythromycin	Methylphenidate	d,l-Propranolol
Aspartame	β-Estradiol	Methylphenidate	d-Propoxyphene
Atropine	Estrorene-3-sulfate	Methylprylon	Pseudoephedrine
Benzic acid	Ethanol (Ethyl alcohol)	Methaqualone	Quinacrine
Benzoic acid	Ethyl-p-aminobenzoate	Metoprolol	Quinidine
Benzoylcegonine	Etodolac	Morphine sulfate	Quinine
Benzphetamine	Famprofazone	Morphine	Ranitidine
Bilirubin	Fenfluramine	3-β-D-glucuronide	Riboflavin
(±)-Brompheniramine	Fenpropfen	Nalidixic acid	Salicylic acid
Bupropion	Fentanyl	Nalorphine	Secobarbital
Caffeine	Fluoxetine	Naloxone	Serotonin
Cannabidiol	Furosemide	Naltrexone	(5-hydroxytryptamine)
Cannabitol	Genisteic acid	Sodium chloride	Sodium thionide
Chloral hydrate	d (+) Glucose	Methypylon	Sulfamethazine
Chloramphenicol	Guaiacol glyceryl ether	Metoprolol	Sulindac
Chlorazepoxide	Guaiacol glyceryl ether carbamate	Nimesulide	Tamazepam
Chloroquine	Hemoglobin	Norcocaine	Tetracycline
Chlorothiazide	Hydralazine	Morphine sulfate	Tetrahydrocortisone, 3-acetate
(+)Chlorpheniramine	Hydrochlorothiazide	α-Naphthaleneacetic acid	Tetrahydrozoline
(±)Chlorpheniramine	Hydrocodone	Norethindrone	Thebaine
Chlorpromazine	Hydrocortisone	Normorphine	Theophylline
Chlorprothixene	Hydromorphone	d-Norpropoxyphene	Thiamine
Cholesterol	p-Hydroxyamphetamine	Orphenadrine	Thioridazine (chlorpromazine)
Cimetidine	o-Hydroxyhippuric acid	Noscapine	I-Thyroxine
Clomipramine	p-Hydroxymethamphetamine	Oxalic acid	Tolbutamide
Clonidine	Hydroxyzine	Oxycodone	cis-Tramadol
Cocaine	lupropfen	Oxymetazoline	Trazodone
Cortisone	Imipramine	Oxymorphone	Triamterene
Creatinine	Iproniazid	Papaverine	Trifluoperazine
Cyclobarbitol	Cyclobenzaprine	Pemoline	Trimethobenzamide
Cyclobenzaprine	Deoxycorticosterone	(-) Deoxyephedrine	Trimethoprim
Dextromethorphan	R (-) Deprenyl	Kanamycin	Trimipramine
Diazepam	R (-) Deprenyl	Pentazocine	Tryptamine
Diclofenac	Diethylpropion	Perphenazine	d,l-Tryptophan
Dicyclomine	Diethylpropion	Phencyclidine	Tyramine
Diffusal	Lidocaine	Phenelzine	d,l-Tyrosine
Dioxin	Lindane	Pheniramine	Uric acid
	(hexachlorocyclohexane)	Phenobarbital	Verapamil
		Phenothiazine	Zomepirac

BIBLIOGRAPHY

- Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*. 6th Edition. Biomedical Publications, Foster City, CA. 2002; 744-747
- Hardman JG, Limbird LE. *Goodman and Gilman's: The Pharmacological Basis for Therapeutics*. 10th Edition. McGraw Hill Medical Publishing, 2001; 208-209.

Index of Symbols

	Attention, see instructions for use		Tests per kit		Authorized Representative
	For in vitro diagnostic use only		Use by		Do not reuse
	Store between 2-30°C		Lot Number		Catalog #
	Do not use if package is damaged				

SCREEN ITALIA S.r.l.
Via dell'Artigianato, 16
06089 - Torgiano - Perugia - Italia
www.screenitalia.it info@screenitalia.it

Número: 145081401
In vigore dal: 2019-11-20